

PHYSICS

CHEMISTRY
BIOLOGY

ENGINEERING



554 801 X-ray apparatus





554 801

X-ray apparatus Mo, complete

Fully-featured, microprocessor-controlled device with x-ray tube Mo and goniometer designed for conducting a wide variety of experiments in x-ray physics. The high-voltage system, x-ray tube and experiment chamber are all within a radiation-proof housing. German type approval as school x-ray apparatus and full-protection device. The type approval is also valid for further x-ray tubes (Fe, Cu, Ag, W). The x-ray tubes are delivered completely adjusted and allow thus an easy and user-friendly exchange. Highest safety and operation comfort by an automatic door locking, which unlock the doors automatically, when no x-ray radiation is generated.

Two large displays show all relevant information on the current experiment. The tube voltage and tube current can be set in the ranges 0 to 35 kV and 0 to 1 mA respectively. The built-in rate meter including counter-tube voltage supply enables direct measuring in conjunction with a Geiger-Müller counter tube. The x-ray apparatus can also be connected to a PC via the USB-port (software included) for recording Bragg spectra. Alternative the two analogue outputs (counting rate and angular position) permit data acquisition using a chart recorder.

The goniometer (554 831) enables precise setting of any angular position of the sensor and the target as well as coupled 2:1 sensor and target motion, both manual and automatic angular scans are possible.

Two screened coaxial lead-ins and one free access duct provide access to set-ups in the experiment chamber, e.g. for connecting an x-ray energy detector.

Device fully assembled and adjusted, ready for operation.

554 800

X-ray apparatus

Basic device fully assembled and adjusted for all tubes, however, without tubes and goniometer.

Scope of delivery: X-ray apparatus without tube, cover for fluorescent screen, dust cover, USB cable, software for Windows 98/2000/XP/Vista

554 980 en

**Book: X-ray apparatus**

24 Experiments from the series "LD Physics Leaflets", english

School x-ray apparatus and full-protection device with German type approval for school use (approval No. BFS 05/07 V/SchRöV) (suitable for the operation with the exchangeable tubes: Fe, Cu, Mo, Ag, W)

Dose rate at a distance of 10 cm: <1 µS/h

Two independent safety circuits for doors, high voltage and emission current

Automatic door locking : doors can be opened only, when no high voltage is present i.e. no x-ray radiation can be generated

High voltage: 0 ... 35.0 kV (regulated DC voltage)

Tube current: 0 ... 1.00 mA (independent regulated DC)

Visible x-ray tube with molybdenum anode for characteristic short-wave radiation:

$K_{\alpha} = 17.4 \text{ keV (71.0 pm)}$, $K_{\beta} = 19.6 \text{ keV (63.1 pm)}$

Fluorescent screen for transillumination experiments: $d = 15 \text{ cm}$

Built-in rate meter including voltage supply for GM counter tube

Loudspeaker: as an acoustic ratemeter

Two 4-digit displays (25 mm high) for displaying the following as desired: high voltage, anode current, counting rate, target/sensor angle, scanning range, step width, gate time

Goniometer (554 831), stepping-motor controlled

Operating modes: manual control and automatic scan for sensor only, target only, 2:1 coupling

Angular range:

target unlimited ($0^{\circ} \dots 360^{\circ}$)

sensor $-10^{\circ} \dots +170^{\circ}$ Step width: 0.1°

Exposure timer, gate time: 0.5 s ... 9999 s

Bushings in the experiment chamber: high-voltage coaxial cable, BNC coaxial cable, empty channel for e.g. tubing, cable etc.

Analog outputs: each proportional to target angle and to counting rate for chart recorder connection

USB port for connecting a PC to control the x-ray apparatus, data recording and evaluation by the delivered Windows software

LabVIEW™ driver for Windows and Linux available free of charge at

<http://www.ld-didactic.com> for user defined controlling and measuring

Input voltage: 230 V ($\pm 10\%$) / 47 - 63 Hz

Power consumption: 120 VA

Dimensions: 67 cm x 48 cm x 35 cm

Weight: 41 kg

Scope of delivery

X-ray apparatus with X-ray tube Mo

Goniometer (554 831)

NaCl crystal (554 78), Lattice-plane spacing: 282 pm

Zirconium foil

Cover for fluorescent screen

Dust cover

USB cable

Software for Windows 98/2000/XP/Vista

Additionally required for 554 800 or 554 801:

End-window counter with cable 559 01



Including Windows software for recording and evaluating data



LabVIEW™

is a registered trademark of National Instruments


X-ray tubes

Directly heated hot cathode tube with screw thread for heat sink and two-pin plug socket for cathode heating for X-ray apparatus (554 800/801)


Max. emission current: 1 mA
Max. anode voltage: 35 kV
Size of focal spot: approx. 2 mm²
Minimum service life: 300 hours
Diameter: 4.5 cm
Length: 20 cm
Weight: 0.3 kg




Bragg-spectrum (NaCl) of a molybdenum tube




Bragg-spectrum (NaCl) of a copper tube



Bragg-spectrum (LiF) of an iron tube with fine structure of the K_{α} -line in 2nd order resolved



Energy spectrum (LiF) of a tungsten tube, in 2nd order 4 L_{β} -lines are well resolved



Bragg-spectrum (NaCl) of a silver tube


X-ray tube	Mo	Cu	Fe	W	Ag
Anode material	Molybdenum	Copper	Iron	Tungsten	Silver
characteristic radiation	$K_{\alpha} = 71.1 \text{ pm}$ (17.4 keV) $K_{\beta} = 63.1 \text{ pm}$ (19.6 keV)	$K_{\alpha} = 154 \text{ pm}$ (8.04 keV) $K_{\beta} = 139 \text{ pm}$ (8.91 keV)	$K_{\alpha} = 194 \text{ pm}$ (6.40 keV) $K_{\beta} = 176 \text{ pm}$ (7.06 keV)	$L_{\alpha} = 148 \text{ pm}$ (8.39 keV) $L_{\beta} = 128 \text{ pm}$ (9.67 keV)	$K_{\alpha} = 56,1 \text{ pm}$ (22.1 keV) $K_{\beta} = 49.7 \text{ pm}$ (24.9 keV)
Absorber foil (to generate monochromatic radiation)	Zirconium (Zr)	Nickel (Ni)	-	-	-


554 831

Goniometer

With two independently controllable stepping motors which move the sensor and target arm. The motion is defined using the keys in the control panel of the X-ray apparatus (554 800 and 554 801) and initiated manually or automatically. Included in the scope of supply of the X-ray apparatus (554 801).

Working principle: stepping motors for target and sensor arms, which can be electronically coupled
Angular range of target: unlimited (0° ... 360°)
Angular range of sensor: approx. -10° to +170°
Angular resolution: 0.1°
Length of sensor arm: approx. 40 -110 mm
Width of sensor slit: 1 mm
Area of target platform: 25 mm x 28 mm
Sample clamping width: 3 - 9 mm
Dimensions: 13.5 cm x 22.5 cm x 12.5 cm
Weight: 3 kg






554 77

LiF crystal for Bragg reflection

Designed to fit the goniometer of the x-ray apparatus (554 801). For experiments in Bragg's configuration, e.g. diffraction (up to the 5th order), x-ray spectra, wavelength determination, Duane and Hunt's displacement law, determining Planck's constant, dependence of absorption on wavelength, determination of lattice plane spacings.

Dimensions: 25 mm x 25 mm x 4 mm
Lattice-plane spacing: 201 pm
Reflection angle for molybdenum K_{α} -radiation (1st order): 10.2°
Crystal structure: face-centered cubic
Surface: parallel [100]



554 842

Set of 2 crystal powder holders

For pressing a crystal powder and then measuring the X-ray diffraction spectrums with powder samples in the X-ray apparatus (554 801).

Dimensions: 25 x 25 x 3 mm each
Weight: 10 g



554 836

Compton accessory X-ray

For x-ray apparatus (554 801) for investigating the Compton effect by means of wavelength-dependent transmission as a function of the placement of the Cu filter in front of or behind the aluminum scattering body; with aluminum scattering body and copper filter in frame.


Aluminum scattering body: 25 x 25 x 4 mm
Copper filter:
Frame: \varnothing 24 mm x 11 mm
Foil: 10 mm x 0.07 mm

554 78

NaCl crystal for Bragg reflection

Design as (554 77)

Lattice-plane spacing: 282 pm
Reflection angle for molybdenum K_{α} -radiation (1st order): 7.24°




554 834

Absorption accessory X-ray

For x-ray apparatus (554 801). Two absorbers for quantitative investigation of the attenuation of x-rays as a function of the thickness and the atomic number of the absorber.

Thickness graduation of aluminum absorber: 0.5/ 1.0/1.5/ 2.0 / 2.5 and 3.0 mm
Material and atomic number for absorbers of constant thickness (0.5 mm):
polystyrene: $Z = 6$, aluminum: $Z = 13$, iron : $Z = 26$, copper : $Z = 29$,
zirconium : $Z = 40$, silver : $Z = 47$
Dimensions of diaphragm: 2.5 x 15 mm
Diaphragm spacing: 5 mm (approx. 10°)
Dimensions : 40 mm x 35 mm x 8 mm each

Scope of delivery: Absorber set I: varying thickness, same material
Absorber set II: varying materials, constant thickness



554 832

Set of absorber foils for X-ray apparatus (554 800/801)

e.g. for experiments on the λ^3 relationship and Moseley's law; foils mounted in frame for attaching to slit diaphragm collimator or counter-tube holder.

Frame: 24 mm \varnothing x 11 mm, Foils: 10 mm \varnothing .

Z	element	thickness
13	Al aluminium	0.5 mm
26	Fe iron	0.5 mm
29	Cu copper	0.07 mm
40	Zr zirconium	0.05 mm
42	Mo molybdenum	0.1 mm
47	Ag silver	0.05 mm
49	In indium	0.3 mm



554 840

Plate capacitor X-ray

For x-ray apparatus (554 800/801) for measuring the ionization current and determining the ion dose rate; electrical connections via 4 mm safety sockets and BNC socket. 3 plug-in mounting pins for defined setup of plate capacitor in the experiment zone.

Input voltage: 0 to 500 V DC
Saturation current: max. 3×10^{-9} A
Saturation voltage: approx. 100 V DC
Ionizable volume of air: 121 cm³
Plate width: 8.5 cm/14 cm
Plate spacing: 3.5 cm
Dimensions: 19 cm x 14 cm x 17 cm

554 839

Blood vessel model for contrast medium

For the demonstration of the effect of contrast media. Plastic plate with covered channels, via screw connections the contrast medium can be injected from outside the x-ray apparatus and its penetration can be observed on the fluorescence screen of the x-ray apparatus. By variation of the distance, magnification effects can be demonstrated.

Scope of delivery:
Plate with blood vessel model on magnet support
Hose
2 plastic syringes
2 stoppers

necessary accessories:
672 6610 potassium iodide, 100 g

554 838


Film holder X-ray

For x-ray apparatus (554 800/801), with printed scale for defined positioning of films for transillumination, Laue diagrams and Debye-Scherrer photographs; includes experiment rail with millimeter scale and pinhole diaphragm $D = 1$ mm for attaching to slit-diaphragm collimator.

Suitable X-ray films:

X-ray film 554 895
X-ray film Agfa Dentus M2 554 896

Dimensions:
Film holder: 12 cm x 16.5 cm
Experiment rail: 25 cm x 16 cm x 6 mm
Pinhole diaphragm: 1 mm Ø



554 8391

Implant model (without picture)

Wood quader with inserted hidden steel pin for the transillumination in the x-ray apparatus.

Additionally required:
Film holder X-ray 554 838



Picture of the blood vessel model (554 839) on fluorescent screen

Crystals for Laue diagrams

554 87 Lithium fluoride crystal
554 88 Sodium chloride crystal

554 895

X-ray film (without picture)

High-sensitivity films for beta, gamma and X-rays, every single sheet is light-sealed and manufactured with developer and fixer. For processing in daylight.

Packet contents: 25 films
Film size: 30 x 40 mm

554 896

X-ray film Agfa Dentus M2 (without picture)

X-ray film welded in light proof plastic foil for use in day light. The film must be removed from the foil for development, for example with the aid of the changing bag (554 899).

Packet contents: 25 films
Film size: 5 cm x 7 cm

554 897

Developer for X-ray film (without picture)

Pack of 10 portions of 125 ml

554 898

Fixative for X-ray film (without picture)

Pack of 10 portions of 125 ml

554 893

Development tank 500 ml (without picture)

For developing the x-ray film (554 896) and for up to two 35 mm films.

554 899

Changing bag for the developer tank (without picture)

Made of double-layer special material. For putting a film in the (554 893) developer tank during daylight hours.

Dimensions: 55 cm x 65 cm

X-ray energy detector

The energy detector is designed for use with the X-ray apparatus (--- 554 801), recording energy separated measurements within the X-ray spectrum. The detector expands the range of experiments possible in quantum and atomic physics, and will also perform non-destructive material analysis.

Topics

Recording the spectrum of an x-ray tube and dependency on current and voltage

Quantitative investigation of the Compton effect

X-ray fluorescence and recording fluorescence spectra of different elements

Verification of Moseley's law using fluorescence spectra

Non-destructive material analysis


Material Analysis:

Geography (rock analysis)

Ecology (detection of heavy elements in nature)

Biology (chemical elements in food)

Chemistry (chemical analysis)



559 938

X-ray energy detector

For the insert in the X-ray apparatus (--- 554 801) for recording of energy dissolved X-ray spectra in connection with Sensor-CASSY (--- 524 010) and MCA box (--- 524 058). The detector contains a thermoelectric cold silicon PIN-detector as well as the electronics for amplification and preparation of the voltage impulses. The amount of the output impuls is proportional with the energy of the X-ray photon.

Photosensitive area: 0.8 mm Ø

Cooling of the detector: thermoelectric (Peltier element)

Entrance window (plastics): absorption equivalent to graphite with $d = 40 \mu\text{m}$

Detectable energy field: approx. 2 keV to 60 keV

Energy resolution at $E = 6.40 \text{ keV}$ (Fe K_{α} -line): 0.4 keV half-width value


Distribution voltage: $\pm 15 \text{ V}$, $+5 \text{ V}$ (via plug-in power supply, in scope of delivery)

Output: BNC socket for connection to the MCA box

Dimensions: 60 mm x 120 mm x 60 mm

Weight: 450 g






Direct determination of the energy loss of the scattered molybdenum K_{α} -radiation (Comptoneffect).

554 8371

Compton accessory X-ray II

For investigating the Compton-effect on X-ray radiation in combination with the X-ray energy detector (559 938) and the X-ray apparatus (554 801). Consists of a circular collimator and a Plexiglas radiation body.

Dimensions: 25 mm x 25 mm x 6 mm





Energy loss of K-lines as a function of the atomic number (Moseley law)

554 844

Set of targets K-line fluorescence

For recording the X-ray fluorescence spectrum of different elements in the X-ray apparatus (554 801) in connection with the X-ray energy detector (559 938), Sensor-CASSY (524 010) and MCA box (524 058).

Materials: Ti, Fe, Ni, Cu, Zn, Zr, Mo, Ag
Dimensions: 25 mm x 25 mm





Excitation of the L-line fluorescence of lead

554 846

Set of targets L-line fluorescence

For recording the X-ray fluorescence spectrum of different elements in the X-ray apparatus (554 801) in connection with the X-ray energy detector (559 938), Sensor-CASSY (524 010) and MCA box (524 058).

Materials: Ag, In, Sn, W, Au, Pb
Dimensions: 25 mm x 25 mm




Determining the chemical composition of alloys by means of x-ray fluorescence analysis


554 848

Set of targets alloys

For recording the X-ray fluorescence spectrum of different alloys in the X-ray apparatus (554 801) in connection with the X-ray energy detector (559 938), Sensor-CASSY (524 010) and MCA box (524 058).

Materials: special steel X5CrNi18-10, brass CuZn36, brass CuZn39Pb3, cobalt samarium magnet
Dimensions: 25 mm x 25 mm





www.ld-didactic.com

Meet the experts: the LD Didactic GmbH Website

Visit our website for the latest special offers and exciting information about our products. Plus, our Online Center is waiting to help you with technical questions.



+ 49 2233 604 - 0

LD Didactic GmbH
Leyboldstraße 1
D – 50354 Huerth - Germany

+ 49 2233 604 - 222

info@ld-didactic.de

www.ld-didactic.com

